



THE NEW DERMAL FILLERS FROM
THE ITALIAN RESEARCH



aurum strong aurum extra strong

Cross Linked Hyaluronic Acid

AURUM fillers are sterile medical devices based on cross-linked and intercalated hyaluronic acid (BDDE), produced by bacterial fermentation. They are injectable, monophasic, apyrogenic, absorbable.

Thanks to the high concentration of HA and a special production process is guaranteed a high performance, a balanced and continuous distribution, an adaptation to the individual morphology. This translates into a full patient satisfaction.



THE NEW DERMAL FILLERS FROM THE ITALIAN RESEARCH


Aurum Strong 24mg

The high elasticity and the medium viscosity of the Aurum strong make it particularly suitable for creating a certain degree of volume and eversion of the lips.

Aurum Extrastrong 26mg


The high concentration and high viscosity of Aurum extrastrong make it particularly suitable for:

- volume recovery
- correction of the cheekbone area
- correction of the nasolabial area
- correction of buccal corners
- correction of the chin area
- reduction of medium and deep wrinkles



TECHNICAL DATA	AURUM STRONG	AURUM EXTRA STRONG
HA concentration	2.4% (24mg)	2.6% (26mg)
Duration	6-8 months	8-10 months
Packaging	1 x 1ml syringe	1x1ml syringe
pH	7	7
Cross-link rate	++++	+++++
BDDE residue	<2PPM	<2PPM
Needle	2 x 27 g	2 x 25 g

Novel Round-Linked Hyaluronic Acid



AURUM fillers are produced with the new ROUND-LINKED HA process. This innovative cross-linking technique allows a less use of BDDE chemical agent guaranteeing long lasting results with greater safety.